

# Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

(This safety data sheet is for information only as it does not comply with the official language requirements of Article 31 (5) of REACH nor does it provide the national information in sections 8 and 15 as specified in Annex II of REACH.)

## dolomite

Version number: 2.0  
Replaces version of: 2023-01-19 (1)

Revision: 2025-04-04  
First version: 2023-01-19

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

<b>Identification of the substance</b>	dolomite
<b>Trade name</b>	<b><u>dolomite</u></b>
<b>Product number</b>	FLX-RAW 01c, FLX-RAW 01b, GQB-01, GQB2-01, CS-0043-CP, FLX-2004
<b>Registration number (REACH)</b>	The substance is exempted from the obligation to register.
<b>EC number</b>	240-440-2
<b>CAS number</b>	16389-88-1

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

<b>Relevant identified uses</b>	Laboratory chemicals
---------------------------------	----------------------

#### 1.3 Details of the supplier of the safety data sheet

FLUXANA® GmbH & Co. KG Borschelstraße 3 D-47551 Bedburg-Hau Germany	Telephone: +49 (0) 2821 - 48011-10 Telefax: +49 (0) 2821 - 48011-99 e-mail: info@fluxana.de Website: www.fluxana.de
--	--

<b>e-mail (competent person)</b>	sdb@csb-compliance.com
----------------------------------	------------------------

Please do not use this e-mail address to ask for the latest safety data sheet. For this purpose contact FLUXANA® GmbH & Co. KG.

#### 1.4 Emergency telephone number

As above or nearest toxicological information centre.

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

##### **Classification according to Regulation (EC) No 1272/2008 (CLP)**

This substance does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

## 2.2 Label elements

### Labelling according to Regulation (EC) No 1272/2008 (CLP)

Not required.

## 2.3 Other hazards

### Results of PBT and vPvB assessment

According to the results of its assessment, this substance is not a PBT or a vPvB.

### Endocrine disrupting properties

Not listed.

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

<b>Name of substance</b>	dolomite
<b>Identifiers</b>	
CAS No	16389-88-1
EC No	240-440-2
<b>Molecular formula</b>	CaMg(CO <sub>3</sub> ) <sub>2</sub>
<b>Molar mass</b>	184.4 g/mol

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

Self-protection of the first aider.

In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.

#### Following skin contact

Wash with plenty of soap and water.

After contact with skin, take off immediately all contaminated clothing, and wash immediately with plenty of water and soap.

#### Following eye contact

Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

#### Following ingestion

Rinse mouth. Do not induce vomiting.

Get medical advice/attention if you feel unwell.

#### Notes for the doctor

None.

## 4.2 Most important symptoms and effects, both acute and delayed

This information is not available.

## 4.3 Indication of any immediate medical attention and special treatment needed

None.

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing media

co-ordinate firefighting measures to the fire surroundings

### 5.2 Special hazards arising from the substance or mixture

Hazardous decomposition products: Section 10.

#### Hazardous combustion products

carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>)

### 5.3 Advice for firefighters

In case of fire and/or explosion do not breathe fumes.

Collect contaminated firefighting water separately.

Fight fire with normal precautions from a reasonable distance.

#### Special protective equipment for firefighters

Self-contained breathing apparatus (EN 133)

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Ventilate affected area.

Do not breathe dust.

Control of dust.

Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

#### For emergency responders

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

### 6.2 Environmental precautions

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

### 6.3 Methods and material for containment and cleaning up

#### Advice on how to contain a spill

Covering of drains.

Take up mechanically.

#### Advice on how to clean up a spill

Take up mechanically.

## **Appropriate containment techniques**

Neutralisation techniques.

## **Other information relating to spills and releases**

Place in appropriate containers for disposal.

Ventilate affected area.

## **6.4 Reference to other sections**

Hazardous combustion products: see section 5.

Personal protective equipment: see section 8.

Incompatible materials: see section 10.

Disposal considerations: see section 13.

## **SECTION 7: Handling and storage**

### **7.1 Precautions for safe handling**

Avoid contact with eyes.

Do not breathe dust.

#### **Measures to prevent fire as well as aerosol and dust generation**

Use local and general ventilation.

Removal of dust deposits.

#### **Handling of incompatible substances or mixtures**

Do not mix with acids.

#### **Measures to protect the environment**

Avoid release to the environment.

Do not empty into drains; dispose of this material and its container at hazardous or special waste collection point.

#### **Advice on general occupational hygiene**

Do not eat, drink and smoke in work areas.

Wash hands after use.

Preventive skin protection (barrier creams/ointments) is recommended.

Remove contaminated clothing and protective equipment before entering eating areas.

### **7.2 Conditions for safe storage, including any incompatibilities**

#### **Flammability hazards**

None.

#### **Incompatible substances or mixtures**

Incompatible materials: see section 10.

#### **Consideration of other advice**

Keep away from food, drink and animal feedingstuffs.

#### **Ventilation requirements**

Provision of sufficient ventilation.

#### **Specific designs for storage rooms or vessels**

Keep container tightly closed and in a well-ventilated place.

**Packaging compatibilities**

Keep only in original container.

**7.3 Specific end use(s)**

No information available.

**SECTION 8: Exposure controls/personal protection****8.1 Control parameters****Occupational exposure limit values (Workplace Exposure Limits)**

This information is not available.

**8.2 Exposure controls****Appropriate engineering controls**

Use local and general ventilation.

**Individual protection measures (personal protective equipment)****Eye/face protection**

Use safety goggle with side protection.

**Hand protection**

Protective gloves		
Material	Material thickness	Breakthrough times of the glove material
these information are not available	these information are not available	these information are not available

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

**Body protection**

Protective clothing for use against solid particulates.

(EN 13832, EN 340, EN 13034, EN 14605).

**Respiratory protection**

In case of inadequate ventilation wear respiratory protection.

Particle filter device (DIN EN 143).

**Environmental exposure controls**

Use appropriate container to avoid environmental contamination.

Keep away from drains, surface and ground water.

**SECTION 9: Physical and chemical properties****9.1 Information on basic physical and chemical properties****Physical state**

solid

**Colour**

grey - white

# dolomite

Version number: 2.0

Revision: 2025-04-04

<b>Odour</b>	odourless
<b>Melting point/freezing point</b>	not determined
<b>Boiling point or initial boiling point and boiling range</b>	not determined
<b>Flammability</b>	non-combustible
<b>Lower and upper explosion limit</b>	not applicable (solid)
<b>Flash point</b>	not applicable
<b>Auto-ignition temperature</b>	not applicable (solid)
<b>Decomposition temperature</b>	>730 °C
<b>pH (value)</b>	9.5 – 10.5 (DIN ISO 787 (7))
<b>Viscosity</b>	not relevant (solid)
<b>Solubility(ies)</b>	
Water solubility	not miscible in any proportion
<b>Partition coefficient n-octanol/water (log value)</b>	not relevant (inorganic)
<b>Vapour pressure</b>	not determined
<b>Density and/or relative density</b>	
Density	2.8 – 2.99 g/cm <sup>3</sup> at 20 °C
Relative vapour density	not relevant (solid)
<b>Particle characteristics</b>	no data available
<b>9.2 Other information</b>	
<b>Information with regard to physical hazard classes</b>	hazard classes acc. to GHS (physical hazards): not relevant
<b>Other safety characteristics</b>	there is no additional information

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

**10.2 Chemical stability**

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

See below "Conditions to avoid".

**10.3 Possibility of hazardous reactions**

No known hazardous reactions.

**10.4 Conditions to avoid**

Control of dust.

**10.5 Incompatible materials**

acids

**10.6 Hazardous decomposition products**

Hazardous combustion products: see section 5.

Carbon dioxide (CO<sub>2</sub>).

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**

If not otherwise specified the classification is based on:

Animal studies; Evidence from any other toxicity tests; Expert judgement (weight of evidence determination).

**Classification according to GHS (1272/2008/EC, CLP)**

This substance does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

**Acute toxicity**

Shall not be classified as acutely toxic (oral).

Exposure route	Endpoint	Value	Species
oral	LD50	6,450 mg/kg	rat

**Skin corrosion/irritation**

Shall not be classified as corrosive/irritant to skin.

**Serious eye damage/eye irritation**

Shall not be classified as seriously damaging to the eye or eye irritant.

**Respiratory or skin sensitisation****Skin sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Respiratory sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Germ cell mutagenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Carcinogenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Reproductive toxicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Specific target organ toxicity - single exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Specific target organ toxicity - repeated exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

**Aspiration hazard**

Shall not be classified as presenting an aspiration hazard.

**11.2 Information on other hazards****Endocrine disrupting properties**

Not listed.

**SECTION 12: Ecological information****12.1 Toxicity****Aquatic toxicity (acute)**

No data available.

**Aquatic toxicity (chronic)**

No data available.

**12.2 Persistence and degradability****Biodegradation**

The study does not need to be conducted because the substance is inorganic.

**Persistence**

The study does not need to be conducted because the substance is inorganic.

**12.3 Bioaccumulative potential****n-octanol/water (log KOW)**

not relevant  
(inorganic)

**12.4 Mobility in soil**

No data available.

## 12.5 Results of PBT and vPvB assessment

According to the results of its assessment, this substance is not a PBT or a vPvB.

## 12.6 Endocrine disrupting properties

Not listed.

## 12.7 Other adverse effects

Data are not available.

### Remarks

Wassergefährdungsklasse, WGK (water hazard class): Nwg.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local/regional/national/international regulations.

#### Sewage disposal-relevant information

Do not empty into drains.

#### Waste treatment of containers/packagings

Completely emptied packages can be recycled.

Handle contaminated packages in the same way as the substance itself.

### Remarks

Please consider the relevant national or regional provisions.

## SECTION 14: Transport information

14.1	UN number or ID number	not subject to transport regulations
14.2	UN proper shipping name	-
14.3	Transport hazard class(es)	-
14.4	Packing group	-
14.5	Environmental hazards	-
14.6	Special precautions for user	-
14.7	Maritime transport in bulk according to IMO instruments	-

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Relevant provisions of the European Union (EU)

#### Restrictions according to REACH, Annex XVII

Not listed.

# dolomite

Version number: 2.0

Revision: 2025-04-04

## List of substances subject to authorisation (REACH, Annex XIV) / SVHC - candidate list

Not listed.

## Seveso Directive

Not assigned.

## Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)

Not listed.

## Regulation on the marketing and use of explosives precursors

Not listed.

## Regulation on drug precursors

Not listed.

## Regulation on substances that deplete the ozone layer (ODS)

Not listed.

## Regulation concerning the export and import of hazardous chemicals (PIC)

Not listed.

## Regulation on persistent organic pollutants (POP)

Not listed.

## 15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out for this substance by the supplier.

## SECTION 16: Other information

### Indication of changes (revised safety data sheet)

Section	Former entry (text/value)	Actual entry (text/value)
1.1	Trade name: GQB-01, RAW01b	Trade name: dolomite
1.1	Product number: FLX-RAW 01b, GQB-01, GQB2-01, CS-0043-CP, FLX-2004	Product number: FLX-RAW 01c, FLX-RAW 01b, GQB-01, GQB2-01, CS-0043-CP, FLX-2004
8.2	-	Body protection: Protective clothing for use against solid particulates. (EN 13832, EN 340, EN 13034, EN 14605).

### Abbreviations and acronyms

Abbr.	Descriptions of used abbreviations
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord relatif au transport international des marchandises dangereuses par route (Agreement con-

# dolomite

Version number: 2.0

Revision: 2025-04-04

Abbr.	Descriptions of used abbreviations
	cerning the International Carriage of Dangerous Goods by Road)
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
EC No	The EC Inventory (EINECS, ELINCS and the NLP-list) is the source for the seven-digit EC number, an identifier of substances commercially available within the EU (European Union)
EINECS	European Inventory of Existing Commercial Chemical Substances
ELINCS	European List of Notified Chemical Substances
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
IMDG	International Maritime Dangerous Goods Code
LD50	Lethal Dose 50 %: the LD50 corresponds to the dose of a tested substance causing 50 % lethality during a specified time interval
NLP	No-Longer Polymer
PBT	Persistent, Bioaccumulative and Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
SVHC	Substance of Very High Concern
vPvB	Very Persistent and very Bioaccumulative

## Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.

Regulation (EC) No. 1907/2006 (REACH), amended by 2020/878/EU.

Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).

International Maritime Dangerous Goods Code (IMDG).

Dangerous Goods Regulations (DGR) for the air transport (IATA).

## Responsible for the safety data sheet

C.S.B. GmbH  
Dujardinstr. 5  
47829 Krefeld  
Germany

Telephone: +49 (0) 2151 - 652086 - 0  
Telefax: +49 (0) 2151 - 652086 - 9  
e-Mail: [info@csb-compliance.com](mailto:info@csb-compliance.com)  
Website: [www.csb-compliance.com](http://www.csb-compliance.com)

## Disclaimer

This information is based upon the present state of our knowledge.

This SDS has been compiled and is solely intended for this product.